

Canadian Aquatic Biomonitoring Network (CABIN) Field Equipment Checklist

General Equipment	
Field sheets and clipboard	
Pencils and markers	
Gloves (rubber, neoprene)	
Waterproof labels	
Labelling tape	
Ziploc bags	
Duct tape and tool kit	
Location and Reach data	
GPS	
Camera	
Densimeter	
Channel and Substrate characteristics	
Velocity metre OR Meter stick	
Measuring Tape	
15 or 30cm ruler	
Hand Level	
Calculator	
Tent pegs	
Water chemistry sampling	
Water quality metres (Temp, pH, DO, Conductivity, turbidity)	
Cooler with sample bottles and ice pack	
Extra batteries	
Benthic Sampling	
Kicknet	
Stopwatch	
Sieve	
White tray	
Squeeze Bottle	
Spoon/tweezers	
Bucket	
Sample jars	
Formalin with MSDS, gloves and glasses	
Cooler for sample jars & Formalin	
Safety equipment	
Lifejackets	
First aid kits (field and vehicle)	
Cell phone or Satellite phone	
Swift water helmet	
Throw bags	
Waders, boots, raingear	
Sunscreen, hat, bug spray	

